



Rust Remover

- Liquid formulation for removing rust
- Inhibits further oxidation
- Suitable for use on iron, steel and all ferrous metals



USAGE Interior	APPLICATION Cloth	DRYING TIME Minimum 24 hours	TOOL CLEANING Water
--------------------------	-----------------------------	---	-------------------------------

Product use

Liberon Rust Remover deoxidises metal surfaces, removing rust and creating a barrier against oxidation. It is suitable for use on ferrous metals, brass, copper and bronze. It can also be used to degrease and clean most metal surfaces prior to colouring with Liberon Antiquing Fluid, Haematite, Iron Paste or paint.

Preparation

Remove loose rust with a wire brush.

Application

Shake before use. Apply Rust Remover with a brush to the affected area only. Allow the surface to whiten, brush off white film with a wire brush and repeat the process until no more white deposit appear. Then apply a final coat to inhibit re-oxidation.

Important: Always test product on a spare surface or inconspicuous area for colour compatibility and end result.

Precautions

- Always wear rubber gloves and eye protection when using this product.

Packaging

Rust Remover is available in 125 ml and 500 ml containers.

Health and Safety

Irritant. Contains Phosphoric Acid. Flammable. Irritating to eyes and skin. Keep out of the reach of children. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not empty into drains. Wear suitable gloves. If swallowed seek medical advice immediately and show this container or label. Use only in well-ventilated areas.